

ASAL
ORIGINAL

PIHAK BERKUASA
PERANTI PERUBATAN



MEDICAL DEVICE
AUTHORITY

PIHAK BERKUASA PERANTI PERUBATAN
MEDICAL DEVICE AUTHORITY
AKTA PERANTI PERUBATAN 2012 (AKTA 737)
MEDICAL DEVICE ACT 2012 (ACT 737)
SIJIL PENDAFTARAN PERANTI PERUBATAN
MEDICAL DEVICE REGISTRATION CERTIFICATE
Seksyen 5(1) Akta 737
Section 5(1) of Act 737

No. Pendaftaran: **GC10845426-228971**
Registration No.:

Tarikh Sah Pendaftaran:
Registration Validity Date:

16/03/2026 - 15/03/2031

Sijil ini adalah dengan ini diberi kepada:
This certificate is hereby issued to:

RAINBOW OCULAR DEVICES SDN BHD

yang beralamat di:
which is located at:

**78A, JALAN BP7/8 BANDAR BUKIT PUCHONG,
47120 PUCHONG, SELANGOR,
47120 PUCHONG
SELANGOR DARUL EHSAN**

Peranan establismen
Role of establishment

Wakil Diberi Kuasa
Authorized Representative

bagi mengesahkan peranti perubatan seperti yang dinyatakan dalam Lampiran 1 adalah berdaftar di bawah Seksyen 5(1) Akta 737.

to confirm that the medical device as detailed out in Attachment 1 is registered under Section 5(1) of Act 737.

Pendaftaran ini diberikan tertakluk kepada peruntukan-peruntukan di bawah Akta 737 dan peraturan-peraturan yang dibuat dibawahnya serta syarat-syarat seperti di Lampiran 2.

This registration is granted subject to the provisions under Act 737 and its subsidiary legislations and the conditions as in Attachment 2.



MURALITHARAN PARAMASUA
KETUA EKSEKUTIF
CHIEF EXECUTIVE
PIHAK BERKUASA PERANTI PERUBATAN
MEDICAL DEVICE AUTHORITY



No. Pendaftaran: **GC10845426-228971**
Registration No.:

Tarikh Sah Pendaftaran:
Registration Validity Date:

16/03/2026 -
15/03/2031

Butir-butir peranti perubatan yang didaftarkan
Particulars of the registered medical device

Nama Peranti Perubatan **PMMA INTRAOCULAR LENSES**
Medical Device Name

Kelas **CLASS C**
Class

Jenama **FIC5585H**
Brand

Kelompok **FAMILY**
Group

Nama dan alamat pembuat:
Name and address of manufacturer

FREEDOM OPHTHALMIC PVT LTD
PLOT NO 31-PHASE-1 SIPCOT INDUSTRIAL COMPLEX, HOSUR -635 126,
TAMIL NADU, INDIA,
635126
INDIA

APPENDIX

NO	NAME AS PER DEVICE LABEL	IDENTIFIER	BRIEF DESCRIPTION OF ITEM
1	PMMA Iris Claw Intraocular Lens	FIC 4272H	Optic Shape, Optic Size & Length
2	PMMA Iris Claw Intraocular Lens	FIC5085H	Optic Shape, Optic Size & Length
3	PMMA Iris Claw Intraocular Lens	FIC5585H	Optic Shape, Optic Size & Length
4	PMMA Posterier Anterior Intraocular Lens	PCC503	Optic Shape, Optic Size & Length
5	PMMA Posterier Anterior Intraocular Lens	PMC524	Optic Shape, Optic Size & Length
6	PMMA Posterier Anterior Intraocular Lens	PMC553	Optic Shape, Optic Size & Length
7	PMMA Posterier Anterior Intraocular Lens	PCC603H	Optic Shape, Optic Size & Length
8	PMMA Posterier Anterior Intraocular Lens	PSF651	Optic Shape, Optic Size & Length
9	PMMA Posterier Anterior Intraocular Lens	PMC 602	Optic Shape, Optic Size & Length
10	PMMA Posterier Anterior Intraocular Lens	PMC602H	Optic Shape, Optic Size & Length
11	PMMA Posterier Anterior Intraocular Lens	PMC652H	Optic Shape, Optic Size & Length



NO	NAME AS PER DEVICE LABEL	IDENTIFIER	BRIEF DESCRIPTION OF ITEM
12	PMMA Posterier Anterior Intraocular Lens	PMS603	Optic Shape, Optic Size & Length
13	PMMA Posterier Anterior Intraocular Lens	PMS602	Optic Shape, Optic Size & Length
14	PMMA 360 Square Adge Aspheric Intraocular Lens	PMC603HSQA	Optic Shape, Optic Size & Length
15	PMMA 360 Square Adge Aspheric Intraocular Lens	PMC603HSQYA	Optic Shape, Optic Size & Length
"End Of Product List"			





MANUFACTURING

📍 Plot # 31, Phase 1, SIPCOT, Hosur - 635 126,
Tamil Nadu, INDIA

☎ + (91)-95855 99420

✉ customercare@freedomophthalmic.com

INTERNATIONAL SALES

📍 134 rue JEAN MERMOZ,
29480 Le Relecq Kerhuon, FRANCE

☎ +33 652940382 📱 +91 95855 99837

✉ fred@freedomophthalmic.com




INTRAOCULAR LENSES

www.freedomophthalmic.com

FREEDOM LENS - POSTERIOR & ANTERIOR CHAMBER



PMMA Posterior & Anterior Chamber Single Piece Intraocular Lens
Lente Intraocular de una Sola Pieza Para Cámara Posterior & Anterior PMMA
Lentille Intraoculaire Monobloc à Chambre Postérieure et Antérieure en PMMA



Model	PCC 503	PMC 524	PMC 553	PCC 603H	PSF 651
Optic Design	Equiconvex				
Optic Size	5.00 mm	5.25 mm	5.50 mm	6.00 mm	6.50 mm
Over All length	12.50 mm	12.00 mm	12.50 mm	12.50 mm	13.50 mm
Haptic Design	Cap. C step vault	Mod. C step vault		Cap. C step vault	Modified C
A-Constant	118.4				118.2
AC Depth*	5.00 mm				4.80 mm




ANTERIOR CHAMBER					
Model	PMC 602	PMC 602H	PMC 652H	PMS 603	PMS 602
Optic Design	Equiconvex			Equiconvex	
Optic Size	6.00 mm		6.50 mm	6.00 mm	
Over All length	13.00 mm			12.50 mm	13.00 mm
Haptic Design	Mod. C step vault			Modified S	
A-Constant	118.2			115.3	
AC Depth*	4.80 mm			2.95 mm	

FREEDOM POWER & GOLD



PMMA 360° Square Edge Aspheric Single Piece Intraocular Lens
Lente Intraocular Plegable Asférico con Borde Cuadrado 360° de PMMA
Lentille Intraoculaire Asphérique Souple en PMMA à Bord Carré sur 360°



Model	PMC 603HSQA	PMC 603HSQYA
Optic Design	Equiconvex - 360° Square Edge with Aspheric Optic	
Optic Size	6.00 mm	6.00 mm
Over All length	12.50 mm	12.50 mm
Haptic Design	Mod. C step vault	
A-Constant	118.4	
AC Depth*	5.00mm	

FREEDOM IRIS

PMMA Iris Claw Intraocular Lens
Lente Intraocular PMMA Iris Claw
Lentille Intraoculaire PMMA Iris Claw



Model	FIC 4272H	FIC 5085H	FIC 5585H
Optic Design	Biconvex - Convex Concave (according to the diopter)		
Optic Size	4.20 mm	5.00 mm	5.50 mm
Over All length	7.20 mm	8.50 mm	8.50 mm
Haptic Design	Claw		
A-Constant	Anterior Side - 115.0 & Posterior Side - 117.4		
AC Depth*	2.95		

FREEDOM RING



PMMA Capsular & Modified Capsular Tension Ring
Anillos de Tensión Capsular & Capsular Modificado PMMA
Anneau de Tension capsulaire & capsulaire modifié en PMMA



Model	FCR 1109	FCR 1210	FCR 1311	FCR 1412
Normal Size	11.00 mm	12.00 mm	13.00 mm	14.00 mm
Compressed Size	9.00 mm	10.00 mm	11.00 mm	12.00 mm



Model	FCC 1210L	FCC 1311L	FCC 1210R	FCC 1311R	FCC 1210LR	FCC 1311LR
Normal Size	12.00 mm	13.00 mm	12.00 mm	13.00 mm	12.00 mm	13.00 mm
Compressed Size	10.00 mm	11.00 mm	10.00 mm	11.00 mm	10.00 mm	11.00 mm

FREEDOM SEGMENT RING



Model	FCS - 10
Length	9.00 mm
Fixation Hook Length	0.40 mm