

OPHTHALMIC CANNULAS



Anaesthesia Cannula Disposable ophthalmic cannula used for administering local anesthesia.



- EO Sterile
- CE certified

Product Code	Description with Gauge	Needle Length	Qty. per box
MIPL/E6/A4	Retrobulbar Cannula (Atkinson) - 19G	38 mm	10 Pcs. (Sterile)
MIPL/E6/A5	Retrobulbar Cannula (Atkinson) - 21G	38 mm	10 Pcs. (Sterile)
MIPL/E6/A1	Retrobulbar Cannula (Atkinson) - 23G	38 mm	10 Pcs. (Sterile)
MIPL/E6/A2	Retrobulbar Cannula (Atkinson) - 25G	38 mm	10 Pcs. (Sterile)
MIPL/E6/A3	Retrobulbar Cannula (Atkinson) - 27G	38 mm	10 Pcs. (Sterile)




Retrobulbar Cannula

Product Code	Description with Gauge	Needle Length	Qty. per box
MIPL/E6/B4	Peribulbar Cannula (Atkinson) - 19G	25 mm	10 Pcs. (Sterile)
MIPL/E6/B5	Peribulbar Cannula (Atkinson) - 21G	25 mm	10 Pcs. (Sterile)
MIPL/E6/B1	Peribulbar Cannula (Atkinson) - 23G	25 mm	10 Pcs. (Sterile)
MIPL/E6/B2	Peribulbar Cannula (Atkinson) - 25G	25 mm	10 Pcs. (Sterile)
MIPL/E6/B3	Peribulbar Cannula (Atkinson) - 27G	25 mm	10 Pcs. (Sterile)



Peribulbar Cannula

Product Code	Description with Gauge	Needle Length	Qty. per box
MIPL/E6/C1	Sub Tenon Cannula - 19G	25 mm	10 Pcs. (Sterile)



Sub Tenon Cannula

OPHTHALMIC CANNULAS




Cystotome Cannula



Disposable ophthalmic cannula used for capsulorhexis, intercapsular / endocapsular and can opener capsulotomy technique. Sharp edge on both sides of the cannula tip is used to puncture and cut/tear the capsule.


- EO Sterile
- CE certified

Product Code	Description with Gauge	Needle Length	Qty. per box
MIPL/E7/A1	Cystotome (Formed) - 25G	16 mm	10 Pcs. (Sterile)
MIPL/E7/A5	Cystotome (Formed) - 26G	16 mm	10 Pcs. (Sterile)
MIPL/E7/A2	Cystotome (Formed) - 27G	16 mm	10 Pcs. (Sterile)
MIPL/E7/A3	Cystotome (Formed) - 30G	16 mm	10 Pcs. (Sterile)



Cystotome (formed)

Product Code	Description with Gauge	Needle Length	Qty. per box
MIPL/E7/A4	Cystotome (Blumenthal) - 27G	16 mm	10 Pcs. (Sterile)




Cystotome (Blumenthal)

Hydrodissection/ Hydrodelineation Cannula

Disposable ophthalmic cannula used to effectively separate layers of nucleus and cortex in cataract surgery procedures.



- EO Sterile
- CE certified



Hydrodissection/
Hydrodelineation Cannula

Product Code	Description with Gauge	Needle Length	Bend Length	Qty. per box
MIPL/E9/A3	Hydrodissection/Hydrodelineation Cannula - 23G	22 mm	9 mm	10 Pcs. (Sterile)
MIPL/E9/A1	Hydrodissection/Hydrodelineation Cannula - 25G	22 mm	9 mm	10 Pcs. (Sterile)
MIPL/E9/A4	Hydrodissection/Hydrodelineation Cannula - 27G	22 mm	9 mm	10 Pcs. (Sterile)
MIPL/E9/A2	Hydrodissection/Hydrodelineation Cannula (J-Shaped) - 25G	22 mm	9 mm	10 Pcs. (Sterile)

OPHTHALMIC CANNULAS



Air Injection & Irrigation Cannula

Disposable ophthalmic cannula used for maintaining or forming anterior chamber and irrigating cornea during surgery. It may be used for injection or removal of air, fluids & viscoelastic substances during an ophthalmic surgery.




- EO Sterile
- CE certified

Product Code	Description with Gauge	Needle Length	Bend Length	Qty. per box
MIPL/E8/A5	Anterior Chamber Cannula (Rycroft) - 16G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A6	Anterior Chamber Cannula (Rycroft) - 19G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A7	Anterior Chamber Cannula (Rycroft) - 20G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A8	Anterior Chamber Cannula (Rycroft) - 21G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A1	Anterior Chamber Cannula (Rycroft) - 23G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A2	Anterior Chamber Cannula (Rycroft) - 25G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A9	Anterior Chamber Cannula (Rycroft) - 26G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A3	Anterior Chamber Cannula (Rycroft) - 27G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A10	Anterior Chamber Cannula (Rycroft) - 28G	22 mm	6 mm	10 Pcs. (Sterile)
MIPL/E8/A4	Anterior Chamber Cannula (Rycroft) - 30G	22 mm	6 mm	10 Pcs. (Sterile)



Anterior Chamber Cannula

Product Code	Description with Gauge	Needle Length	Bend Length	Qty. per box
MIPL/E8/B1	Viscoelastic Cannula - 23G	22 mm	9 mm	10 Pcs. (Sterile)
MIPL/E8/B2	Viscoelastic Cannula - 25G	22 mm	9 mm	10 Pcs. (Sterile)
MIPL/E8/B3	Viscoelastic Cannula - 27G	22 mm	9 mm	10 Pcs. (Sterile)




Viscoelastic Cannula

Lasik Cannula

Disposable ophthalmic cannula used for all types of refractive surgery and minor office procedure.

Product Code	Description with Gauge	Needle Length	Qty. per box
MIPL/E16/A1	Lasik Irrigating Cannula (Banaji) - 23G	22 mm	10 Pcs. (Sterile)
MIPL/E16/A2	Lasik Double Arm Cannula - Irrigation by eight ports	N/A	10 Pcs. (Sterile)



Lasik Cannula

- EO Sterile


OPHTHALMIC CANNULAS



Cortex Aspiring Cannula

Disposable ophthalmic cannula used for aspirating cortical debris when chamber is maintained with either an anterior chamber maintainer or a viscoelastic substance.

Product Code	Description with Gauge	Needle Length	Qty. per box
MIPL/E10/A1	Cortex Aspiring Cannula (Simcoe) - 21G	22 mm	10 Pcs. (Sterile)
MIPL/E10/A2	Cortex Aspiring Cannula (Simcoe) - 23G	22 mm	10 Pcs. (Sterile)
MIPL/E10/A4	Cortex Aspiring Cannula (Simcoe) - 25G	22 mm	10 Pcs. (Sterile)
MIPL/E10/A5	Cortex Aspiring Cannula (Simcoe) - 26G	22 mm	10 Pcs. (Sterile)
MIPL/E10/A3	Cortex Aspiring Cannula (J-Shaped) - 25G	22 mm	10 Pcs. (Sterile)



Cortex Aspiring Cannula

- EO Sterile
- CE certified

Irrigating & Aspiring Cannula

Disposable ophthalmic cannula used for maintaining or forming anterior chamber and irrigating cornea during surgery. It may be used for injection or removal of air, fluids & viscoelastic substances during an ophthalmic surgery.

Product Code	Description with Gauge	CE Certified	Qty. per box
MIPL/E11/A6	I & A Cannula (Simcoe Regular) - 21G	No	10 Pcs. (Sterile)
MIPL/E11/A1	I & A Cannula (Simcoe Regular) - 22G	Yes	10 Pcs. (Sterile)
MIPL/E11/A2	I & A Cannula (Simcoe Regular) - 23G	Yes	10 Pcs. (Sterile)



Irrigating & Aspiring Cannula

- EO Sterile

Capsule Polishers

Disposable ophthalmic cannula used for removing debris from the posterior, anterior and equatorial regions of capsular bag. Micro etched tips are slightly roughened for gentle polishing.

Product Code	Description with Gauge	Needle Length	Bend Length	Qty. per box
MIPL/E12/A1	Capsule Polisher (Kratz) - 21G	22 mm	7 mm	10 Pcs. (Sterile)
MIPL/E12/A2	Capsule Polisher (Kratz) - 23G	22 mm	7 mm	10 Pcs. (Sterile)
MIPL/E12/A4	Capsule Polisher (Kratz) - 25G	22 mm	7 mm	10 Pcs. (Sterile)
MIPL/E12/A5	Capsule Polisher (Kratz) - 27G	22 mm	7 mm	10 Pcs. (Sterile)
MIPL/E12/A3	Capsule Polisher (Pearce) Olive Tip - 25G	25 mm	7 mm	10 Pcs. (Sterile)



Capsule Polishers

- EO Sterile
- CE certified

OPHTHALMIC CANNULAS



Bimanual Cannula

Disposable ophthalmic cannula used for irrigation and aspiration of the anterior chamber.



- EO Sterile
- CE certified

Bimanual I & A Cannula - Disposable

Product Code	Description	Gauge	Port Size	Qty. per box
MIPL/E15/DP/21WP	Bimanual Cannula Handpiece Set with Polisher – Disposable			1 Set (Sterile)
MIPL/E15/A8/01	Bimanual Irrigation Handpiece – Disposable	21G	0.50mm	1 Pc. (Sterile)
MIPL/E15/A8/02	Bimanual Aspiration Handpiece with polisher – Disposable	21G	0.40mm	1 Pc. (Sterile)
MIPL/E15/DP/22WP	Bimanual Cannula Handpiece Set with Polisher – Disposable			1 Set (Sterile)
MIPL/E15/A9/01	Bimanual Irrigation Handpiece – Disposable	22G	0.45mm	1 Pc. (Sterile)
MIPL/E15/A9/02	Bimanual Aspiration Handpiece with polisher – Disposable	22G	0.35mm	1 Pc. (Sterile)
MIPL/E15/DP/23WP	Bimanual Cannula Handpiece Set with Polisher – Disposable			1 Set (Sterile)
MIPL/E15/A10/01	Bimanual Irrigation Handpiece – Disposable	23G	0.40mm	1 Pc. (Sterile)
MIPL/E15/A10/02	Bimanual Aspiration Handpiece with polisher – Disposable	23G	0.30mm	1 Pc. (Sterile)



- Supplied Non-Sterile

Bimanual I & A Cannula - Titanium

Product Code	Description	Gauge	Port Size	Qty. per pack
MIPL/E15/A6/01	Bimanual Irrigation Handpiece - Titanium	21G	0.50mm	1 Pc.
MIPL/E15/A6/02	Bimanual Aspiration Handpiece with polisher - Titanium	21G	0.40mm	1 Pc.
MIPL/E15/A7/01	Bimanual Irrigation Handpiece - Titanium	22G	0.45mm	1 Pc.
MIPL/E15/A7/02	Bimanual Aspiration Handpiece with polisher - Titanium	22G	0.35mm	1 Pc.
MIPL/E15/A5/01	Bimanual Irrigation Handpiece - Titanium	23G	0.40mm	1 Pc.
MIPL/E15/A5/02	Bimanual Aspiration Handpiece with polisher - Titanium	23G	0.30mm	1 Pc.

Coaxial Cannula

Disposable ophthalmic cannula used for irrigation and aspiration of the anterior chamber.



- EO Sterile

Coaxial I & A Handpiece - Disposable

Product Code	Description	Gauge	Port Size	Qty. per pack
MIPL/E19/A2/01	Coaxial I&A Handpiece (Straight) - Disposable	21G	0.40mm	1 pc. (Sterile)
MIPL/E19/A2/02	Coaxial I&A Handpiece (45° Angle) - Disposable	21G	0.40mm	1 pc. (Sterile)
MIPL/E19/A3/01	Coaxial I&A Handpiece (Straight) - Disposable	23G	0.30mm	1 pc. (Sterile)
MIPL/E19/A3/02	Coaxial I&A Handpiece (45° Angle) - Disposable	23G	0.30mm	1 pc. (Sterile)