

Rainbow Ocular Devices Sdn Bhd

Co No: 202401048133 (1593978-A)

DISTRIBUTOR'S OF
OPHTHALMOLOGY PRODUCTS

Bringing Clearer
Vision to Life



Precision
Diagnostics



Trusted
Products



Better
Patient Care

COMPANY PROFILE 2026



Reg. No 222 25003



Reg. No 81-306-5442



Reg. No MDA-8226-W125



Malaysia Ministry of Finance (MOF)
Registered Company

Reg. No 357-000-2422982



+60-102405768



03-8066-7448, 010-266 1852



www.rainbowoculardevices.com



info@rainbowoculardevices.com



ABOUT US

Rainbow Ocular Devices Sdn Bhd was established in November 2024 in Malaysia. It is dedicated to supplying high-quality and affordable ophthalmology products throughout Malaysia. Our primary mission is to support eye care professionals and institutions by providing reliable, innovative, and customer-focused solutions. Over time, we have earned our reputation as a trusted partner within the ophthalmic sector by maintaining our commitment to quality, continuous innovation, and customer satisfaction.

The company was founded by Dr. Thanga Babu, an experienced ophthalmologist, which uniquely positions Rainbow Ocular Devices as both a potential end user and a vendor. This dual perspective enables us to rigorously assess the products we offer. If a product meets our quality checks, we are confident that it will meet the high standards of the global ophthalmology community as well.

At this moment I would like to thank Dato Paduka Zarina Sabdin, our esteemed CORPORATE ADVISOR for being in our board and guiding us.

Dr.Thanga Babu MD (Eye) MBA (IB)
Managing Director

Dato. Paduka Zarina Sabdin
Corporate Advisor

Ms. Archana Thanga Babu
M.Sc. Biomedical Sciences, Clinical Nutritionist
CEO The Rainbow Group

PRODUCT'S PORTFOLIO

Our products are certified by the Malaysian Medical Device Authority (MDA), ensuring compliance with stringent regulatory standards.


Diagnostic Strips

We represent Madhu Instruments India, offering a comprehensive range of diagnostic strips, with a focus on Fluorescein Strips and Schirmer's Strips and its Blue Variant.



OVD's (Viscoelastics)

We distribute Dispersive and Cohesive OVD's. Our partnership with Siloam Ophthalmics, India, allows us to supply their complete range of viscoelastics, including HPMC and Sodium Hyaluronate 1.4.



Trypan Blue

We provide Trypan Blue from Siloam Ophthalmics, India. Concentration of 0.80mgm/mL makes the anterior capsule staining enhanced to do the rhexis procedure hassle free.



Sodium Hyaluronate (3% w/v) with Sodium Chondroitin Sulfate (4% w/v) Ophthalmic Solution – iOcuCoat

iOcuCoat - Ophthalmic Viscosurgical Device has an intermediate cohesive/dispersive index (CDI). It is a sterile, non-pyrogenic, viscoelastic solution of a highly purified non-inflammatory medium molecular weight fraction of sodium hyaluronate and sodium chondroitin sulphate. iOcuCoat solution maintains a deep chamber during anterior segment surgeries, enhances visualization during the surgical procedure, and protects the corneal endothelium and other ocular tissues.



The viscoelasticity of the solution maintains the normal position of the vitreous face, thus preventing formation of a postoperative flat chamber. It is designed to create and maintain space, to protect the corneal endothelium and other intraocular tissues and to manipulate tissues during surgery.

Molecular Weight : 2.8 - 3.2 Million Dalton

Viscosity at Rest : 300000 - 450000 mPas

pH : 6.8 - 7.6

Osmolality: 320— 370 mOsm/Kg

Raw Material: Sodium Hyaluronate - Non-Animal Origin

Sodium Chondroitin Sulfate - Non-Bovine Origin

Advantage: Imported EP Grade API

No refrigeration required - Store at Room Temp

BD Glass Syringe & Aseptic Filling

**Packing Material : Premium Quality of Pharmaceutical Grade.
Available in 1ml PFS**

Eye Shield

Madhu Instruments India supplies us with high-quality, sterile, transparent eye shields. Most economical and wide distributed all over Malaysia



Visco-Adaptive OVD

We represent ODC Italy, for their Co-Visc 2% Visco Adaptive Sodium Hyaluronate OVD contains medium and short chain low molecular weight. This combines the advantages of the cohesive physical properties of the viscoelastic inactive and the dispersive physical properties of low molecular weight.



Ophthalmic Knives

Our collaboration with Shah Ophthalmic Pvt Ltd enables us to offer Acucut Knives, a globally recognised and accredited range of ophthalmic knives. Their full range of knives are used by many ophthalmologists all over Malaysia.

1.0mm to 1.6mm



Easy Knife Identification

- Color-coded handles
- Printed blade tip names

Anti Glare (Matt Finish)


- Improves view of the Operating field
- Eliminate reflective glare from operating room light source
- Maximum visibility of the knife within the eye chamber

Enhanced Control

- Ergonomic design
- Textured non-slip Handle

Depth Indicator

- For Accurate Incisions



Intra-Ocular Lenses (IOLs)

We represent the Optima, India brand in Malaysia, distributing a full suite of IOLs, including Hydrophilic Phakic Toric and Hydrophobic Trifocal models.




ERGOTOUCH FLY BLISTER



- 1x Dispenser box contains 10x blister packs
- 1x Blister pack contains 1x injector and 1x cartridge
- Dispenser box: 8.9cm width, 22.2cm height and 20.4cm depth
- Blister packs: 8.6cm width, 2.3cm height and 20.2cm depth
- Blister packs and dispenser boxes provide extra safety for both the injector and the cartridge, due to the rigid and strong packaging

Eye Drapes Surgical Set


We provide high-quality ophthalmic eye drapes sourced from Biotex, Turkey. Including Eye Surgery Drape Kit which includes eye drape with pouch, drape, back table drape, Standard surgeon gown, Towel and Pack of Sterilization.



	Includings / İçerik	Size / Ebat	Piece / Adet
1	Eye Drape With Pouch / Tek Poşlu Göz Örtüsü	100 * 130 cm	1
2	Drape / Düz Örtü	150 *180 cm	1
3	Back Table Drape / Alet Masa Örtüsü	100 *150 cm	1
4	Standard Surgery Gown / Standart Cerrahi Önlük	Large	2
5	Towel / Havlu	40 * 40 cm	2
6	Pack Of Sterilization / Sterilizasyon Bohçası	100 * 100 cm	1

Eye Drape With Pouch

Product Name	Dimensions	Piece	Fabric Type
Sterile Ophthalmic Drape Single Pouch / Rod-338	40*50 CM	1	40gr SMS
Sterile Ophthalmic Drape Single Pouch / Rod-338-1	100*130 CM	1	40gr SMS
Sterile Ophthalmic Drape Single Pouch / Rod-338-2	100*150 CM	1	40gr SMS
Sterile Ophthalmic Drape Single Pouch / Rod-338-3	130*130 CM	1	40gr SMS
Sterile Ophthalmic Drape Single Pouch / Rod-338-4	150*160 CM	1	40gr SMS
Sterile Ophthalmic Drape Single Pouch / Rod-338-5	120*120CM	1	40gr SMS
Sterile Ophthalmic Drape Single Pouch / Rod-338-6	60*60 CM	1	40gr SMS
Sterile Ophthalmic Drape Single Pouch / Rod-338-7	50*50CM	1	40gr SMS



Iris-Claw Intra Ocular Lenses

Our Iris-Claw IOLs are sourced from Freedom Ophthalmics India, whose products are highly regarded by Malaysian ophthalmologists.




Slit Lamp Biomicroscope

We distribute the complete range of Labomed Slit Lamps from the USA, including advanced imaging slit lamp models. Mag's available are 3 step and 5 step.



Operation Theatre Table

As distributors of Plantech OT Tables from India, we have installed these tables in multiple locations across Malaysia. This chair type OT table has gained attention of many eye centre's due to its ergonomic design and cost efficiency.



Capsular Tension Ring for Cataract Surgery

We distribute Excellent Hi Care capsular tension rings Malaysia to support ophthalmic surgeons managing cataract cases with suspected or confirmed zonular weakness. Available in both manual and preloaded variants, these CTRs are designed to help maintain capsular bag contour, support posterior capsule stability, and assist with intraocular lens centration during phacoemulsification and lens implantation.



Clinical Relevance

In cataract surgery, zonular compromise may increase the risk of capsular instability, posterior capsule rupture, vitreous loss, intraocular lens decentration, and postoperative lens displacement. A capsular tension ring provides circumferential support within the capsular bag, redistributing tension across the equator and helping to stabilize areas affected by weak or deficient zonules.

Product Availability

Excellent Hi Care CTRs are offered in manual and preloaded configurations to support surgeon preference and operating room workflow. The manual variant provides familiar controlled insertion, while the preloaded option is intended to simplify handling and reduce preparation steps before implantation.



Preloaded Optima Vue-Safe

The Preloaded Optima Vue-Safe has gained increasing clinical relevance in Malaysia due to its convenience and growing use among cataract surgeons. By presenting the ring in a ready-to-use delivery format, the system supports efficient intraoperative handling while preserving the core purpose of CTR implantation: maintaining capsular bag stability and supporting predictable IOL positioning in appropriate cases.



PMMA Capsular Tension Ring and Advance Capsular Tension Ring


Model	OCR1109	OCR1210	OCR1311	OCR1412	OCR1210L	OCR1311L
Overall Diameter	11.0 mm	12.0 mm	13.0 mm	14.0 mm	12.0 mm	13.0 mm
Insertion Diameter	09.0 mm	10.0 mm	11.0 mm	12.0 mm	10.0 mm	11.0 mm
Fixation Hook	Nil				Left side with Step Vault	
Material	CQ- PMMA Clear / Blue					



Model	OCR 1210R	OCR 1311R	OCR 1210LR	OCR 1311LR
Overall Diameter	12.0 mm	13.0 mm	12.0 mm	13.0 mm
Insertion Diameter	10.0 mm	11.0 mm	10.0 mm	11.0 mm
Fixation Hook	Right Side with Step Vault		Both Side with Step Vault	
Material	CQ- PMMA Clear / Blue			

Ophthalmic Cannula

We are proud distributors of PRICON, India, supplying their full range of ophthalmic cannula globally accepted brand and most of the OVD manufacturers worldwide use their visco cannula only.




Microsurgical Instruments

From high quality stainless steel instruments to high end Titanium ophthalmic microsurgical instruments, we represent PRICON, India as their Malaysian counterpart.



Sterilization Trays

To protect valuable surgical instruments, we provide Sterixx plastic sterilization trays in all sizes, equipped with silicon mats. Our trays are compatible with all standard sterilisation methods, including steam, dry heat, ETO, and chemical processes.



Iris Retractors

We represent Madhu Instruments India for iris retractors, which are widely distributed and utilised across Malaysia. Our partnership with Madhu Instruments ensures that Malaysian ophthalmologists have access to high-quality iris retractors designed to facilitate optimal surgical outcomes.





Lens Glide

Its disposable and used for insertion of IOL's or other microsurgical instruments in the anterior chamber and it's a double ended polypropylene device where the narrow end is used during phacoemulsification procedure and broader end is used during non- phaco surgical procedures.




Pupil Expansion Redefined

For all mid dilate pupil during cataract surgery we are distributing B-Hex pupil expander all models apart from manual iris retractors from Madhu Instruments. B-Hex and B-OCTA are FDA certified pupil expanders from Med Invent Devices, India.



B-HEX® 23 G Forceps



This instrument is a forceps featuring short, curved shafts and serrated jaws. The jaws open in a plane perpendicular to the squeeze handle, facilitating effective grasping and manipulation of B-HEX and B-OCTA pupil expanders. It is sterilized using moist heat through autoclaving. Available in Stainless steel also.

Ai Enabled Online Perimetry

Melbourne Rapid Fields (MRF) AI software is a highly regarded, portable, and cloud-based visual field (perimetry) testing tool designed as an alternative or complement to traditional, bulky, and expensive in-clinic machines.

Developed by the University of Melbourne and the Royal Victorian Eye and Ear Hospital, it is designed for rapid diagnosis and monitoring of glaucoma, age-related macular degeneration (AMD), and diabetic retinopathy. Our company is proud to be its Malaysian Representative.

We represent MRF Ai Enabled Online Perimeter also for Indonesia, Cambodia, Singapore, Hongkong & philippines.

GLANCE Optical Pty Ltd.

Melbourne Rapid Fields (MRF) AI-enabled online perimeter

100,000+ MRF visual fields
performed globally



TRANSFORMATIVE ONLINE FIELD TESTING AT YOUR FINGERTIPS

Validated

Validated against gold standards (proven accuracy and liability with over 20+ clinical trials performed on 1500+ patients)
Registered with Therapeutic Goods Administration (TGA, ARTG: 282166) Australia, MedSafe (NZ), CDSCO, MHRA (UK) and CE Marked as a perimetry device.



Web Cam

Patented Web Camera Technology - To ensure optimal testing condition (including viewing distance monitoring, background lighting condition, gaze monitoring etc.)



Small Foot Print

low startup and running costs



Device Independent

Web browser based perimetry

Progression Analysis

Unique trend and pointwise progression analysis to monitor and detect glaucoma change



Unique Vision Tests

Designed for detecting visual acuity change in different vision conditions (including low contrast and noise VA)



Secure

Secure HIPPA compliant Data Storage

Clinic and Home Based Use



Comprehensive

Standard testing protocols (24-2, 10-2, 30-2, Binocular Esterman equivalent, 24-2 rapid test, Neurological field test and visual acuity tests).
Voice guidance in 12 languages.
Offline mode for testing without network connection.

In collaboration with



CONTACT US



info@rainbowoculardevices.com
www.rainbowoculardevices.com

VISCOT Specialty Skin Markers

We represent Viscot Medical LLC, USA who provides specialized, sterile, and non-sterile, ultra-fine tip skin markers designed for high-precision ophthalmic and plastic procedures. These markers, such as the 1425SR, use scrub-resistant, non-toxic gentian violet ink to create durable, highly visible, and fine, precise lines on the skin.

Cataract- Toric Lens Marking



1424 TWIN TIP BARREL
 Blephmarker: 2 Ultra Fine Tips |
 Traditional Gentian Violet Ink

PART #	DESCRIPTION	CASE QTY
1424-100	Sterile 2 Ultra Fine Tips	100
1424SR-10	Sterile 2 Ultra Fine Tips w/ Ruler	10
1424SR-100	Sterile 2 Ultra Fine Tips w/ Ruler	100

1436 VALUE SURGICAL MARKER
 Ultra Fine Tip | Value Barrel



PART #	DESCRIPTION	CASE QTY
1436-100	Non Sterile Ultra Fine Tip	100
1436SR-100	Sterile Ultra Fine Tip w/ Ruler	10
1436SR-50	Sterile Ultra Fine Tip w/ Ruler	50
1436SR-10	Sterile Ultra Fine Tip w/ Ruler	10
1436SRL9-100	Sterile Ultra Fine Tip w/ Ruler & 9 Blank Label	100



1402/1403 PRECISION SURGICAL MARKER
 Ultra Fine Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1402-50	Sterile Ultra Fine Tip w/ Ruler	50
1402-100	Sterile Ultra Fine Tip w/ Ruler	100
1403-100	Non-Sterile Ultra Fine Tip	100

1436 VALUE SURGICAL MARKER
 Ultra Fine Tip | Value Barrel



PART #	DESCRIPTION	CASE QTY
1425-100	Sterile 1 Regular & 1 Ultra Fine Tip	100
1425SR-100	Sterile 1 Regular & 1 Ultra Fine Tip w/ Ruler	100

Ophthalmology Diagnostic Device Tables & Surgeon Stool's

We represent Plantech, India, presents a premium selection of ophthalmology equipment tables specially designed for high-value diagnostic instruments. Understanding that your sophisticated ophthalmology devices deserve an equally refined base, our motorised, silent tables are engineered to enhance both the performance and aesthetics of your equipment.

Our comprehensive range includes tables with tabletops of various sizes, ensuring that imaging systems along with their CPUs remain stable and secure. With Plantech's solutions, you can showcase your expensive diagnostic tools on tables that not only complement their value but also offer reliable functionality for daily clinical use.



DIAGNOSTIC & SURGICAL LENSES

Single Mirror Goniolens




This gonio lens is used for viewing the anterior chamber angle and central retina.

Due to its small size, a single mirror goniolens is very convenient to handle. The mirror angle is 62° which when rotated through 360° gives a clear view of the anterior chamber angle. Its small contact surface makes it very useful while diagnosing adults as well as children with small palpebral fissure.

- Made in glass
- Contact medium required

Product Code	Description	Mirror Angle	Qty. per box
MIPL/I1	Single Mirror Goniolens	62°	1 pc.



Two Mirror Goniolens




This gonio lens is used for viewing the anterior chamber angle and central retina.

Two opposing mirrors inclined at 62° provide a complete 360° view of angle with only 90° of rotation of the lens. Similar to a single mirror goniolens, the small contact surface helps in examination of children and adults with small fissures.

- Made in glass
- Contact medium required

Product Code	Description	Mirror Angle	Qty. per box
MIPL/I2	Two Mirror Goniolens	62°/62°	1 pc.



Three Mirror Goniolens




This gonio lens is used for viewing the anterior chamber angle and central & peripheral retina.

This lens has three mirrors angled at 59°, 66° and 73° and arranged at 120° intervals to permit viewing of the fundus and anterior chamber angle. The posterior pole is also observable through the central axis of the lens.

- Made in glass
- Contact medium required

Product Code	Description	Mirror Angle	Qty. per box
MIPL/I3	Three Mirror Goniolens	59°/66°/73°	1 pc.



DIAGNOSTIC & SURGICAL LENSES

Four Mirror
Goniolens


This is a standard gonio lens which is used for static and dynamic gonioscopy.



Four mirrors angled at 62° help in viewing all quadrants of the anterior chamber without rotation. The small contact surface of 9 mm diameter makes the lens very convenient to use. It has a unique detachable handle. Coupling fluid not required.

- Made in glass
- No contact medium required

Product Code	Description	Mirror Angle	Qty. per box
MIPL/I4	Four Mirror Goniolens	62°	1 pc.



20D
Aspheric Lens

This lens is used for general examination of the fundus using “Binocular Indirect Ophthalmoscope(BIO)”.



It provides a high resolution image of the retina in the OPD or the operating room. A 50° field of view helps in visualization up to the mid peripheral region. A dynamic examination with patients eye movements helps in visualization of the peripheral retina as well. A working distance of 50mm makes lens manipulation very easy for the user.

- Made in glass.
- Aspheric Optics.

Product Code	Description	Field of View	Qty. per box
MIPL/K1	20D Aspheric Lens	50°	1 pc.

FASTEST FULLY AUTOMATED FUNDUS CAMERA

- **Imaging in 15 Seconds per eye**
- A fully automated, non-mydratic fundus camera that captures high-quality images with minimal training.


Compact Design

A space-saving system that integrates easily into any setting, bringing advanced eye screening anywhere.



Laser-Based Autofocusing System


Delivers crystal-clear retinal images using AI-powered pupil tracking.



Minimum pupil size : 2.6mm

Delivers high-quality images even through pupil as small as 2.6mm.

CARESORB® POLYGLACTIN 910



PRODUCT INFORMATION



COMPOSITION

90% glycolide
10% L-Lactide



COATING

Polyglactin 370
Calcium stearate



STRUCTURE – COLOR

Multifilament
Violet – Undyed



WOUND SUPPORT

~ 30 days



TOTAL ABSORPTION

56-70 days



RETENTION STRENGTH

75% @ 2 weeks
50% @ 3 weeks



USP SIZES

8/0 – 2



KNOT-PULL TENSILE STRENGTH

Exceed USP
requirements (USP)



SUTURE LENGTH

30 cm, 45 cm
75 cm, 90 cm



STERILISATION

EO
(Ethylene Oxide)



PACKAGING

Aluminum Foil
Box 12 pcs

CARESORB Suture Specification Table

No.	CPT code CARESORB	Competitor Code VICRYL NOVOSYN POLYSORB	Suture Type	USP Size	Suture Length (cm)	Needle Shape	Needle Length (mm)	Circle	How Supplied
1	GT02H06L10	W1702	Polyglactin 910	10/0	10	Spatula	6	3/8	Box / 12 pcs
2	GT02GG06L30	V448G	Polyglactin 910	10/0	30	Spatula x2	6	1/2	Box / 12 pcs
3	GT03H06L15	W1703	Polyglactin 910	9/0	15	Spatula	6	3/8	Box / 12 pcs
4	GT04HH06L30	W9559	Polyglactin 910	8/0	30	Spatula x2	6	3/8	Box / 12 pcs
5	GT04HH06L45	W9560	Polyglactin 910	8/0	45	Spatula x2	6	3/8	Box / 12 pcs
6	GT04GG07L30	W9564	Polyglactin 910	8/0	30	Spatula x2	7	1/2	Box / 12 pcs
7	GT05HH06L30	W9561 L1796K	Polyglactin 910	7/0	30	Spatula x2	6	3/8	Box / 12 pcs
8	GT05GG07L45	W9565	Polyglactin 910	7/0	45	Spatula x2	7	1/2	Box / 12 pcs
9	GT07I08L45	W9552 J555G	Polyglactin 910	6/0	45	Spatula x2	8	1/4	Box / 12 pcs
10	GTU07I08L45	W9752	Polyglactin 910 (U)	6/0	45	Spatula x2	8	1/4	Box / 12 pcs
11	GTU07E12PL45	W9500T Safil C1049607	Polyglactin 910 (U)	6/0	45	Reverse Cutting (P)	12	3/8	Box / 12 pcs
12	GT07EE12PL45	W9537	Polyglactin 910	6/0	45	Reverse Cutting (P) x2	12	3/8	Box / 12 pcs
13	GT10I08L45	W9553	Polyglactin 910	5/0	45	Spatula x2	8	1/4	Box / 36 pcs
14	GTU10I08L45	W9753	Polyglactin 910 (U)	5/0	45	Spatula x2	8	1/4	Box / 36 pcs
15	GTU10E12PL45	W9501T	Polyglactin 910 (U)	5/0	45	Reverse Cutting (P)	12	3/8	Box / 36 pcs
16	GTU10HH12L45	W9761	Polyglactin 910 (U)	5/0	45	Spatula x2	12	3/8	Box / 36 pcs

CARELON® NYLON – POLYAMIDE

PRODUCT INFORMATION

**COMPOSITION**

Polyamide 6/6.6

**COATING**

None

**STRUCTURE – COLOR**Monofilament
Blue – Black**USP SIZES**

10/0 – 2

**KNOT-PULL TENSILE STRENGTH**Exceed USP
requirements (USP)**SUTURE LENGTH**30 cm, 45 cm
75 cm, 100 cm**STERILISATION**EO
(Ethylene Oxide)**PACKAGING**Tyvek Pouch
Box 12 | 24 pcs

CARELON® NYLON – SUTURE SPECIFICATION

No.	CPT Code CARELON	Competitor Code ETHILON DAFILON MONOSOF	Suture Type	USP Size	Suture Length (cm)	Needle Shape	Needle Length (mm)	Circle	How Supplied
1	M01B03L15	W2881 G1118366	Nylon	11/0	15	Taper Point	3	3/8	Box / 12 pcs
2	M02B03L15	G1117041	Nylon	10/0	15	Taper Point	3	3/8	Box / 12 pcs
3	M02B05L15	W2814 W2850 W2810 G1118110	Nylon	10/0	15	Taper Point	5	3/8	Box / 12 pcs
4	M02HH06L30	W1770 7770G G1118706 Alcon 8065711701	Nylon	10/0	30	Spatula x2	6	3/8	Box / 12 pcs
5	M02GG07L30	W1719	Nylon	10/0	30	Spatula x2	7	1/2	Box / 12 pcs
6	M03B05L15	W2829 W2813 G1111434	Nylon	9/0	15	Taper Point	5	3/8	Box / 12 pcs
7	M03HH06L30	W1769 W1717 G1118714 Alcon 8065192201	Nylon	9/0	30	Spatula x2	6	3/8	Box / 12 pcs
8	M04B06L15	W2808 G1118170	Nylon	8/0	15	Taper Point	6	3/8	Box / 12 pcs
9	M04E06L30		Nylon	8/0	30	Reverse Cutting	6	3/8	Box / 12 pcs
10	M05HH08L30		Nylon	7/0	30	Spatula x2	8	3/8	Box / 12 pcs
11	M05E13		Nylon	7/0	75	Reverse Cutting	13	3/8	Box / 12 pcs
12	M05EE13L45		Nylon	7/0	45	Reverse Cutting x2	13	3/8	Box / 12 pcs
13	M07E12PL45	W1610T 697G C0932060 SN3697	Nylon	6/0	45	Reverse Cutting (P)	12	3/8	Box / 12 pcs
14	M07D13		Nylon	6/0	75	Reverse Cutting	13	1/2	Box / 12 pcs
15	M07E13	C0932060 C0934061	Nylon	6/0	75	Reverse Cutting	13	3/8	Box / 12 pcs
16	M07E16	C0935115 SN660	Nylon	6/0	75	Reverse Cutting	16	3/8	Box / 12 pcs
17	M10II08L45	Alcon 8065211501	Nylon	5/0	45	Spatula x2	8	1/4	Box / 24 pcs
18	M10E12PL45	W1611T W1852T	Nylon	5/0	45	Reverse Cutting (P)	12	3/8	Box / 24 pcs
19	M10E13	C0932078 Dermalon 1741-21 8886174121	Nylon	5/0	75	Reverse Cutting	13	3/8	Box / 24 pcs
20	M10E16	C0935123 SN-659 SN659	Nylon	5/0	75	Reverse Cutting	16	3/8	Box / 24 pcs
21	M15E16	W3318 C0935131	Nylon	4/0	75	Reverse Cutting	16	3/8	Box / 24 pcs

SOFT PL SILICONE PUNCTUM PLUGS

DETAILS

- Pen-Style Inserter with dilator tip - 1pk, 2pks, 6pks
- Original Inserter with dilator tip - 1pk, 2pks, 6pks, 20pks
- Sterile, preloaded on inserter in 1pk, 2pks, 6pks dispenser box
- Low-profile cap intended to reduce dislocation
- Economy pack includes 20 plugs. Holders equipped with an alignment slot for quick loading



Original Inserter




Pen-Style Inserter

Original Inserter Ref #	Pen-Style Inserter Ref #	QTY/Box	Product Name	Diameter
6610-S †	7710-S †	1pk	SOFT PLUG Silicone, Preloaded, Micro	0.4mm
6611-S †	7711-S †	1pk	SOFT PLUG Silicone, Preloaded, Mini	0.5mm
6612-S †	7712-S †	1pk	SOFT PLUG Silicone, Preloaded, Petite	0.6mm
6613-S †	7713-S †	1pk	SOFT PLUG Silicone, Preloaded, Small	0.7mm
6614-S †	7714-S †	1pk	SOFT PLUG Silicone, Preloaded, Medium	0.8mm
6610 †	7710 †	2pk	SOFT PLUG Silicone, Preloaded, Micro	0.4mm
6611 †	7711 †	2pk	SOFT PLUG Silicone, Preloaded, Mini	0.5mm
6612 †	7712 †	2pk	SOFT PLUG Silicone, Preloaded, Petite	0.6mm
6613 †	7713 †	2pk	SOFT PLUG Silicone, Preloaded, Small	0.7mm
6614 †	7714 †	2pk	SOFT PLUG Silicone, Preloaded, Medium	0.8mm
6610-D †	7710-D †	6pk	SOFT PLUG Silicone, Preloaded, Micro	0.4mm
6611-D †	7711-D †	6pk	SOFT PLUG Silicone, Preloaded, Mini	0.5mm
6612-D †	7712-D †	6pk	SOFT PLUG Silicone, Preloaded, Petite	0.6mm
6613-D †	7713-D †	6pk	SOFT PLUG Silicone, Preloaded, Small	0.7mm
6614-D †	7714-D †	6pk	SOFT PLUG Silicone, Preloaded, Medium	0.8mm
6610-EP		20pk	SOFT PLUG Silicone, Economy Pack, Micro	0.4mm
6611-EP		20pk	SOFT PLUG Silicone, Economy Pack, Mini	0.5mm
6612-EP		20pk	SOFT PLUG Silicone, Economy Pack, Petite	0.6mm
6613-EP		20pk	SOFT PLUG Silicone, Economy Pack, Small	0.7mm
6614-EP		20pk	SOFT PLUG Silicone, Economy Pack, Medium	0.8mm

SOFT PLUG®
SILICONE PUNCTUM PLUG,
FLOW CONTROL

LONG TERM OCCLUSION WITH PARTIAL DRAINAGE

The distinctive low-profile cap combined with a softer, flexible silicone maximizes patient comfort. The inserter minimizes potential for accidental release of the plug. Plug is designed with a center channel allowing for partial tear flow. SOFT PLUG® Flow Control Plugs are removable with forceps (per instruction for use).



DETAILS


- Available in 1pk and 2pks
 - Silicone
- Sterile, preloaded on inserter
 - For partial occlusion with limited tear drainage

Original Inserter Ref #	Pen Inserter Ref #	QTY/Box	Product Name	Diameter
6612-SFC †		1pk	SOFT PLUG Flow Control, Preloaded, Petite	0.6mm
6613-SFC †		1pk	SOFT PLUG Flow Control, Preloaded, Small	0.7mm
6614-SFC †		1pk	SOFT PLUG Flow Control, Preloaded, Medium	0.8mm
6612-FC †	7712-FC ‡	2pk	SOFT PLUG Flow Control, Preloaded, Petite	0.6mm
6613-FC †	7713-FC ‡	2pk	SOFT PLUG Flow Control, Preloaded, Small	0.7mm
6614-FC †	7714-FC ‡	2pk	SOFT PLUG Flow Control, Preloaded, Medium	0.8mm

SOFT PLUG®
COLLAGEN
INTRACANALICULAR PLUGS

TEMPORARY OCCLUSION

The Collagen Plugs have a smooth proprietary finish that simplifies insertion of the plug. Collagen Plugs are ideal for short-term, diagnostic, and postsurgical occlusion.



DETAILS

- Absorbable material lasts 2-5 days
 - Proprietary finish simplifies insertion
- Evaluate patient acceptance
 - Insert with forceps
- 2mm length
 - Sterile

Ref #	QTY/Box	Product Name	Diameter
6502	60 per box	SOFT PLUG Collagen Intracanalicular Plug	0.2mm
6503	60 per box	SOFT PLUG Collagen Intracanalicular Plug	0.3mm
6504	60 per box	SOFT PLUG Collagen Intracanalicular Plug	0.4mm



Quality Assurance

We source our products exclusively from reputable manufacturers and ensure that each item undergoes rigorous quality control processes. Adhering to both national and international standards, we guarantee the safety, reliability, and efficacy of every product we supply.



Industry Partnerships

Through collaborations with leading global brands and innovative technology partners, we bring the latest advancements in ophthalmology to the Malaysian market. These partnerships enable us to deliver comprehensive solutions tailored to the ever-evolving needs of eye care professionals.



Customer Support

Our dedicated customer service team is always available to provide technical support, detailed product information, and prompt after-sales service. We also offer training programmes, workshops, and educational resources to help our clients stay updated with the latest techniques and products in ophthalmology. Our sales team has over a decade of experience in ophthalmology sales across Malaysia, ensuring knowledgeable and reliable support.



Why Choose Us

Why choose Rainbow Ocular Devices? We offer unmatched product quality and reliability, with a wide range of ophthalmology products under one roof. Our efficient logistics ensure prompt delivery across Malaysia, backed by an experienced team dedicated to customer satisfaction. With the unique advantage of being both an end user and a vendor led by an Ophthalmologist at the helm—we stand apart in the industry.

Enhancing VISION

CONTACT US

OPHTHALMOLOGY PRODUCTS / MEDICAL DEVICES

LEADERSHIP

Dr. Thanga Babu MD(Eye) MBA (IB)

Founder & Managing Director

Dato Paduka Zarina Sabdin

Corporate Advisor

LET'S TALK

T 03-8066-7448 / 010-266 1852

W +60-102405768

@ info@rainbowoculardevices.com

WWW www.rainbowoculardevices.com

A 78 A, Jalan BP 7/8, Bandar Bukit Puchong, 47120 Puchong,
Selangor, MALAYSIA

SCAN TO CONNECT

Choose the channel that works best for you



WHATSAPP



LINKEDIN



RAINBOWOCULARDEVICES

INSTAGRAM