

Sterile Ophthalmic Drapes

With collection pouch

FEATURED SIZES

60 x 60 cm

REF SO-6060

1 sterile drape per pack

120 x 120 cm

REF SO-120120

1 sterile drape per pack



Oval configuration shown for illustration.

A defined surgical configuration

25 x 30 cm

Collection pouch

10 x 10 cm

Oval fenestration

15 x 15 cm

3M incise film

45 GSM SMMS NONWOVEN FABRIC

Spunbond + Meltblown + Meltblown + Spunbond

Shelf life: 5 years from date of manufacture.

Material & product range

SMMS construction

Spunbond
Meltblown
Meltblown
Spunbond

4 layers. 45 GSM.

Nonwoven fabric construction.

Layer diagram is schematic;
thickness is not to scale.

Documented size / code range

Current manufacturing focus: 60 x 60 cm and 120 x 120 cm.

DRAPE SIZE	PRODUCT CODE	FOCUS
60 x 60 cm	SO-6060	Priority size
120 x 120 cm	SO-120120	Priority size
100 x 130 cm	SO-100130	-
100 x 150 cm	SO-100150	-
150 x 160 cm	SO-150160	-

Handling & storage

Single use only. Do not resterilize. Consult instructions for use.

Do not use if package is damaged. Keep dry and away from sunlight.

Keep away from open flame, intense heat and flammable anaesthesia.

MANUFACTURED FOR

SILOAM OPHTHALMICS

Level 1, No:23, Periyar Street, Kannigapuram,
Avadi, Chennai 600054, INDIA

GST: 33AWWPS8520K1ZR